# plasma care®

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REDUCTION OF BACTERIAL LOAD

SAFE – HEALTHY TISSUE IS NOT DAMAGED

SUITABLE FOR PATIENTS WITH PACEMAKERS

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REDUCTION OF  
BACTERIAL LOAD

SAFE –  
HEALTHY TISSUE  
IS NOT DAMAGED

SUITABLE FOR  
PATIENTS WITH  
PACEMAKERS

PORTABLE  
AND  
EASY-TO-USE

ACTIVATION  
OF WOUND  
HEALING

# plasma care®

## case reports

The **plasma care®** is used for patients with chronic wounds of various genesis, e. g. postoperative wound healing disorders, diabetic foot syndrome and decubitus. Here you will find some case reports:

### EXAMPLE: POST-OPERATIVE WOUND HEALING DISORDER

40-year-old male patient. Resection of a liposarcoma and partial resection of the costal arches C2-C5 right and C2-C4 left. Post-operative wound healing disorder with sternal putrid (purulent) wound abscess cavity – wound revision and application of a wound sealing system unsuccessful, formation of necrosis.

#### Treatment course:

- Cold plasma treatment over 13 weeks with an average of 3 applications/week.
- Regular removal of crusts and necroses, necrotic areas became smaller, the wound became flatter and closed from inside and from its edges.

Initial situation, Oct 2019

End of January 2020

### EXAMPLE CASE: DECUBITUS ON THE NECK

64-year-old female patient. Post-polio syndrome, chronic respiratory and ventilatory insufficiency, Gr. III obesity, hypercapnic coma. Hypergranulating pressure ulcer on the neck, due to friction and moisture on the strap of the respiratory mask, stagnation in the granulation phase.

#### Treatment course:

- 7 cold plasma treatments (6 x 2 min, 1 x 1 min) over 3 weeks, first treatments were performed on 6 days in a row.
- The wound was completely epithelialized after 22 days.

(Picture documentation following page)

Continued picture documentation decubitus

Day 0, Initial situation

Day 1, 1 CAP treatment

Day 8, 6 CAP treatments

Day 15, 7 CAP treatments

Day 22, completely epithelialized

### EXAMPLE: POST-OPERATIVE WOUND HEALING DISORDER

77-year-old female patient, secondary healing wound (split skin transplantation and flap wound plastic) as a result of trauma (contusion) in a car accident, wound healing disorder with known CVI and heart failure, infection with enterobacter aerogenes, additional tendency to edema formation in the lower legs.

#### Treatment course:

- 2 cold plasma treatments per week (1 min each).
- changeover to moist wound care to soften the incrustations and necroses
- mechanical cleaning
- after 18 treatments over 11 weeks, wound completely epithelialized except for punctual superficial skin gap (wound size: 0.22 x 0.17 cm).

Week 0, Initial situation

Week 1, 3 CAP Treatments

Week 2, 5 CAP treatments

Week 8, 7 CAP treatments

Week 11, 18 CAP treatments

ONLY 1 MINUTE  
PER APPLICATION

USER-FRIENDLY

NO DEVELOPMENT  
OF RESISTANCES

REDUCTION OF  
BACTERIAL LOAD

SAFE –  
HEALTHY TISSUE  
IS NOT DAMAGED

SUITABLE FOR  
PATIENTS WITH  
PACEMAKERS

PORTABLE  
AND  
EASY-TO-USE

ACTIVATION  
OF WOUND  
HEALING

# plasma care® – small medical device with big effect


The **plasma care®** is on market since July 2019. This medical device uses cold atmospheric plasma to inactivate microorganisms including multi-resistant pathogens. Even stimulation of wound healing can also be observed in some patients. Here are some case reports on wound treatment:

## EXAMPLE: POST-OPERATIVE WOUND HEALING DISORDER

63-year-old female patient with a traumatic transtibial amputation of the right lower leg, post-operative wound healing disorder on the stump (support surface to the prosthesis), followed by a flap plastic surgery, no wound healing for over 1½ years – no infection.

### Treatment course:

- 9 cold plasma treatments (1 min/13 cm<sup>2</sup>) within 3.5 weeks.
- Reduction of wound size from 2.04 cm x 1.87 cm (L x W) to 0.48 cm x 0.51 cm (L x W)

Day 0, Initial situation      Day 1, 1 CAP treatment      Day 7, 3 CAP treatments      Day 21, 6 CAP treatments      Day 25, 9 CAP treatments

## EXAMPLE CASE: DIABETIC FOOT SYNDROME

50-year-old male patient with diabetic foot syndrome. Infected ulcers above the metatarsal bones (left) and necrotic 4th toe, amputation due to worsening of the wound situation, wound infection with partially resistant corynebacteria, enterococci and staphylococci.

### Treatment course:

- 2 cold plasma treatments per week for the first 3 weeks, followed by application every 14 days thereafter (9 treatments in 12 weeks) as part of the dressing change
- Healing within 12 weeks after beginning of plasma therapy

(Picture documentation following page)

Continued picture documentation diabetic foot syndrome


Day 0, Initial situation      Day 2, 2 CAP treatment      Day 14, 4 CAP treatments      Day 42, 7 CAP treatments      Day 83, 9 CAP treatments

## EXCAMPLE CASE: ULCUS CRURIS

77-year-old female patient, bedridden due to a cervical spinal injury. Ulcus recurrence of unclear origin on the right lower leg, no edema, vascular structure unclear. Wound healing stagnated for several months, partly purulent coatings, patient suffers from severe pain during mechanical wound cleaning.

### Treatment course:

- 8 cold plasma treatments (1 min) within 4 weeks resulted in significant reduction of wound size
- Progressive epithelialization, pain relief
- Interruption of CAP treatment caused ulcer recurrence
- Complete epithalization of wound 4 weeks after restart of CAP therapy (2 x 1 min/week)

Initial situation      8 CAP treatments      Relapse after interruption of therapy      restart: after 8 CAP treatments

ONLY 1 MINUTE  
PER APPLICATION

USER-FRIENDLY

NO DEVELOPMENT  
OF RESISTANCES

REDUCTION OF  
BACTERIAL LOAD

SAFE –  
HEALTHY TISSUE  
IS NOT DAMAGED

SUITABLE FOR  
PATIENTS WITH  
PACEMAKERS

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REDUCTION OF  
BACTERIAL LOAD

SAFE –  
HEALTHY TISSUE  
IS NOT DAMAGED

SUITABLE FOR  
PATIENTS WITH  
PACEMAKERS

PORTABLE  
AND  
EASY-TO-USE

ACTIVATION  
OF WOUND  
HEALING

# Notes

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

terraplasma  
MEDICAL

**terraplasma medical GmbH**

Parking 32 | 85748 Garching (Munich)  
Germany

Phone.: + 49 89 588 055 30

Fax: + 49 89 588 055 99

info@terraplasma-medical.com

www.terraplasma-medical.com

**Jens Kirsch** | Chief Executive Officer (CEO)

**Petra Speidel** | Chief Financial Officer (CFO)

**PD Dr. rer. nat. habil. med. Julia Zimmermann** | Chief Medical Officer (CMO)

**Prof. Dr. Dr. h.c. Gregor Morfill** | Scientific advisory council

ONLY 1 MINUTE  
PER APPLICATION

USER-FRIENDLY

NO DEVELOPMENT  
OF RESISTANCES

REDUCTION OF  
BACTERIAL LOAD